

PROMOTION VALID 4/01/19 THROUGH 6/30/19

BUY TWO



BUY TWO SETS OF  
CETYLCIDE G

(#0122)

AND RECEIVE ONE  
HYDE-OUT

(#0128)

ABSOLUTELY FREE



GET  
ONE  
FREE!

## CETYLCIDE G

Easy to measure  
convenient tip-and-pour bottle

No waste  
activate only as much as you need

Simple  
formulate with ordinary tap water

Cetylcode-G is a glutaraldehyde-based sterilant/high-level disinfectant. Each quart of concentrate yields 15 quarts (3.75 gallons) of 3.2% strength use-solution, that may be reused for up to 28 days. An 8-oz. bottle of diluent is included with each quart of Cetylcode-G.

## HYDE-OUT NEUTRALIZER

Fast-acting  
5 minutes to complete neutralization

No-mess dispensing  
easy tip-and-pour bottle

Effective  
reduces harmful vapors, too

Hyde-Out neutralizes glutaraldehyde and OPA high level disinfectant/sterilant solutions and converts the solutions into non-toxic, non-hazardous waste. Each bottle neutralizes one gallon of glutaraldehyde or eight gallons of Ortho-phthalaldehyde (OPA). Complies with CA Health and Safety Code section 25123.5 (c).

Free product offers are available from the manufacturer after purchasing qualifying products through your dental supplies distributor. Valid only in the USA. Each offer is final and may not be combined with any other offers. Limit one of each offer per customer/practice. Promotion valid 4/01/19 through 6/30/19. Cetylite must receive redemption request by 7/15/19.

To redeem free product, forward a copy of your original paid dealer invoice, showing purchase of qualifying products within promotion dates.

Email: [redemptions@cetylite.com](mailto:redemptions@cetylite.com)

Mail: Cetylite, Inc., Attn: Customer Service  
9051 River Road, Pennsauken, NJ 08110

Fax: 856-665-5408

Please allow 4-6 weeks for delivery of free goods. Free products to be shipped by Cetylite Industries. Cetylite reserves the right to void any non pre-approved promotion.

9051 River Road  
Pennsauken NJ 08110  
[www.cetylite.com](http://www.cetylite.com)



CETYLITE®