

# 4<sup>TH</sup> QUARTER DENTAL SPECIALS 2018 (USA)

Offers valid October 1, 2018 through December 31, 2018

## SECURE FIT™ MASK TECHNOLOGY

Better fit. Better filtration.  
Better protection.



Aluminum nose and chin strips reduce gapping and fogging.

**ASTM LEVEL 3** High Fluid Protection

50/box

**GCFXSSF-Q418** Ultra™ Sensitive Earloop Mask w/ Secure Fit™ Mask Technology - White

40/box

**GCFXSSFSF-Q418** Ultra™ Sensitive FogFree™ Earloop Mask w/ Secure Fit™ Mask Technology - White

25/box

**GCPWSSF-Q418** Ultra™ Sensitive FogFree™ Earloop Mask w/ Shield w/ Secure Fit™ Mask Technology - White

**ASTM LEVEL 2** Moderate Fluid Protection

50/box

**GCPBLSF-Q418** Procedural Earloop w/ Secure Fit™ Mask Technology - Blue

**GCPVLSF-Q418** Procedural Earloop w/ Secure Fit™ Mask Technology - Lavender

**GCPPKSF-Q418** Procedural Earloop w/ Secure Fit™ Mask Technology - Pink

Learn more: [Crosstex.com/SecureFit](http://Crosstex.com/SecureFit)

Not made with natural rubber latex

**BUY 3  
GET 1**  
shipped with order

## Barrier Film with Finger-Lift Edge



Non-adhesive finger-lift edge allows quick, easy placement/removal. Ideal for areas that are hard to reach, clean or disinfect.

1200 sheets/roll (4" x 6" / 10 x 15 cm)

**BFBL-Q418** Blue  
**BFCL-Q418** Clear  
**BFPP-Q418** Pink with a purpose

Not made with natural rubber latex

**BUY 3  
of the same  
GET 1**  
shipped with order

## Le Petit Masque

Smaller size mask, with fluid resistant outer layer with soft spunbound inner layer.



**ASTM LEVEL 1** Low Fluid Protection

50/box

**GCAPTPK-Q418** Pink

Learn More: [Crosstex.com/SecureFit](http://Crosstex.com/SecureFit)

Not made with natural rubber latex

**BUY 3  
GET 1**  
shipped with order

## Type 4 Sterilization Pouches

Built-in external/internal multiparameter indicators.

200/box

**SCXX2-Q418** 2.25" x 4"  
**SCX2-Q418** 2.75" x 9"  
**SCXS2-Q418** 3.5" x 5.25"  
**SCS2-Q418** 3.5" x 9"  
**SCW2-Q418** 5.25" x 6.5"  
**SCM2-Q418** 5.25" x 10"  
**SCL5152-Q418** 5.25" x 15"  
**SCL2-Q418** 7.5" x 13"  
**SCL8162-Q418** 8" x 16"

100/box

**SCL3222-Q418** 3.5" x 22"  
**SCL10152-Q418** 10" x 15"  
**SCL10172-Q418** 10.5" x 17"  
**SCL12152-Q418** 12" x 15"  
**SCL12182-Q418** 12" x 18"

500/box

**SCS25-Q418** 3.5" x 9"

Learn more: [Crosstex.com/Sure-Check](http://Crosstex.com/Sure-Check)

**BUY 5  
GET 1**  
equal or lesser value item shipped with order



The new look of nitrous oxide/oxygen sedation  
Try NEW Axess Low Profile Nasal Mask

Introductory kit contains:

- 6 single-use disposable sample masks (2 small Bubblegum, 2 medium Mint and 2 large Unscented)
- 1 reusable, sterilizable scavenging circuit, and adapter fittings to retrofit your existing systems.



**52006-Q218**  
Axess™ Scavenger Circuit for RFS™ Remote Flow System and Accu-Vac™ Mount



**52005-Q218**  
Axess™ Scavenger Circuit for in-line controller

**SPECIAL**  
Introductory Offer

## Sure-Check™ Strips

Multiparameter indicator for use in steam and chemical vapor sterilization processes. Use in pouches, tubing, autoclave bags and CSR wrapped trays/containers.

100/bag

**SCK-Q418** White 3.75" x .625" / 10 x 2 cm




**BUY 3  
GET 1**  
shipped with order

# 4<sup>TH</sup> QUARTER DENTAL SPECIALS 2018 (USA)

Offers valid October 1, 2018 through December 31, 2018

## Comfort Plus™ Premium Saliva Ejectors

The kinder, gentler saliva ejector... takes in fluid, not tissue.



### 100/bag

ZWWCP-Q418	White/White
ZCWCP-Q418	Clear/White
ZWPBPG-Q418	White/Pink (Bubblegum Scented)
ZWGCPM-Q418	White/Green (Mint Scented)

Comfort Plus® is a registered trademark of Crosstex International, Inc., a Cantel Medical Company.

**BUY 4**  
of the same  
**GET 1**  
shipped with order

## STEAMPlus™ Type 4 Integrators



Immediate-readout; confirm sterilization conditions with a definitive pass or fail result.

### 100/pack

SSI-100-Q418	4" x .75" Strip
--------------	-----------------

Learn More: [Crosstex.com/SteamPlus](http://Crosstex.com/SteamPlus)

SteamPlus™ is a trademark of SPSmedical Supply Corp.

**BUY 2**  
**GET 1**  
shipped with order

## Sparkle™ Propyl Paste

### Propyl Paste

Advanced formulation improves paste consistency. Spreads easily, rinses completely.



### Coarse Grit 12oz/340gr jar

UPC12-Q418	Mint
------------	------

### Coarse Grit 200/box

UPCB-Q418	Berrylicious Asst.
UPCBG-Q418	Bubble Gum
UPCCH-Q418	Cherry
UPCM-Q418	Mint
UPCOV-Q418	Orange Vanilla
UPCA-Q418	Assorted Flavors

### Medium Grit 200/box

UPMB-Q418	Berrylicious Asst.
UPMBG-Q418	Bubble Gum
UPMCH-Q418	Cherry
UPMM-Q418	Mint
UPMOV-Q418	Orange Vanilla
UPMA-Q418	Assorted Flavors

### Fine Grit 200/box

UPFBG-Q418	Bubble Gum
------------	------------

Assorted flavors include: Cherry, Mint, Orange Vanilla  
Berrylicious assorted flavors include: Raspberry, Blueberry, Strawberry

**BUY 2**  
of the same  
**GET 1**  
shipped with order

## SANI Vac™

Evacuation System Cleaner



A concentrated enzymatic cleaner for continuous, high-performance suction.

### 1-gallon bottle

JVAC-Q418	SANI Vac™ Evacuation System Cleaner
-----------	-------------------------------------

**BUY 1**  
**GET 1**  
shipped with order

## twist™ plus

Less splatter. Less friction.

The only disposable prophyl angle that oscillates in a 120° counter-rotational motion to reduce splatter and frictional heat.



### 100/box

TPLUSPASC-Q418	Soft Cup, gray
TPLUSPAFC-Q418	Firm Cup, white
TPLUSPATB-Q418	Tapered Brush

▶ Watch the video: [Crosstex.com/twistPlus](http://Crosstex.com/twistPlus)

Not made with natural rubber latex. Manufactured in China.

**BUY 2**  
of the same  
**GET 1**  
shipped with order

## Sani-Tab™ Chain-Free™ Towels

The only chain-free patient towel with patented adhesive tabs eliminating the need for bib chains and clips.



### 400/Case

WEXBLT-Q418	Blue
WEXLVT-Q418	Lavender
WEXDRT-Q418	Dusty Rose
WEXWHT-Q418	White

**BUY 4**  
of the same  
**GET 1**  
shipped with order

Buyer's price may constitute a "discount or other reduction in price" under 42 U.S.C. § 1320a-7b(b) and 42 C.F.R. 1001.952(h). Buyer is solely responsible for (i) accurately reporting all discounts and price reductions claimed or charged to any Federal or State health care program, and (ii) determining whether to report or pass on savings or discounts to payors.

### Free goods fulfillment

All offers to be shipped via participating authorized Crosstex dealers, with automatically filled FREE goods. Dealer OTD sales recap required to receive Crosstex free goods credit. Must be sold during valid promotional dates.

Mix & Match Free Goods (must be equal or less value than purchased goods)

Dealer recap reports must be received at Crosstex within 30 days of promotional expiration dates.

### Program details

Promotions MUST be advertised to participate in the Crosstex program (print proof submission required).<sup>1</sup> Dealers may mix and match additional promotions to advertise. Quarterly promotional offers may not be combined with any other offers or promotions.

<sup>1</sup>Proof of purchase invoice