


Dia-ProT™ Gutta Percha

Dia-ProT™ Gutta Percha and Paper Point perfectly and precisely fit the shapes created by Dia-X™ Rotary Files and ProTaper® Rotary Files.

- Millimeter markings provide a quick and easy method for depth measurements
- Cadmium free, latex free, and highly radiopaque
- Length = 28 mm +/- 1
- Convenient one-hand slide open packaging

Buy 8 boxes of Dia-ProT Gutta Percha Points, Get a FREE box of Dia-X™ NiTi Rotary Files of your choice!



Dia-Pro T Gutta Percha (60/box)

- ML150-S601 Dia-ProT F1
- ML150-S602 Dia-ProT F2
- ML150-S603 Dia-ProT F3
- ML150-S604 Dia-ProT F4
- ML150-S605 Dia-ProT F5
- ML150-S691 Dia-ProT Assorted (F1/F2/F3)
- ML150-S692 Dia-ProT Assorted (F4/F5)
- ML150-S693 Dia-ProT Assorted (F1-F5)

{ COMPATIBILITY CHART }

- Dia-X DX: Access File (19mm)
- Dia-X D1: ProTaper® S1
- Dia-X D2: ProTaper® S2

- Dia-X **D3** : ProTaper® F1: Dia-ProT **F1**
- Dia-X **D4** : ProTaper® F2: Dia-ProT **F2**
- Dia-X **D5** : ProTaper® F3: Dia-ProT **F3**

NEW!

Dia-ProT™ NEXT Gutta Percha & Paper Point



Dia-ProT™ NEXT Gutta Percha and Paper Point perfectly and precisely fit the shapes created by ProTaper® NEXT Rotary Files.

- Millimeter markings provide a quick and easy method for depth measurements
- Cadmium free, latex free, and highly radiopaque
- Length = 28 mm +/- 1
- Convenient one-hand slide open packaging

LAUNCH SPECIAL!

DiaProT Next Gutta Percha (60/box)

- ML167-S602 X2
- ML167-S603 X3
- ML167-S604 X4
- ML167-S605 X5
- ML167-S691 Ass't (X2/X5)
- ML167-S692 Ass't (X2/X3)

DiaProT Next Paper Points (100/box)

- MP267-S602 X2
- MP267-S603 X3
- MP267-S604 X4
- MP267-S605 X5
- MP267-S691 Ass't (X2/X5)
- MP267-S692 Ass't (X2/X3)

15% OFF*

Stainless Steel Hand Files



**10
27 OR +1
+ 3***

- Exceeds all ISO/ADA specifications
- Stainless steel files with endo stops
- Color coded to save time
- 6 files/box. Eligible series: 501, 502, 503, 504
- K-Files: available in 21, 25, and 31 mm. Sizes #8-80 and assorted.
- H-Files: available in 21, 25, and 31 mm. Sizes #8-80 and assorted.
- Reamers: available in 21 and 25 mm. Sizes #8-80 and assorted.
- Barbed Broaches: available in 21 and 25 mm. Sizes #20-40 and assorted.



#4007-2000 Dia-Duo Regular Kit

Contains: Duo-Gun Handpiece, Duo-Pen Handpiece & Duo-Charger, Heating Condensers .04-50 & .06-50, Disposable Sheath (200/box), GP Obturation Tips (23G, 4/box), Safety Caps (2/pack), Plunger Set, Plunger O-Ring Assembly (2/pack), Multi Tool, Scrubbing Brush, Power Adapter & Cord

Duo-Gun Cordless Backfill Obturation Device

- Precise extrusion of gutta percha with
- Strong and flexible tips available in 23G and 25G
- Temperature reaches 200°C within 15 seconds
- Ergonomically designed for both left and right-handed individuals with power button and display on both sides
- Display shows battery life and temperature setting
- Variable temperature settings (140°C, 170°C, 200°C)
- Operate hand trigger with minimal effort
- Fast full charge in 90 min. Operation time of 15 x 5 min.

DIA-DUO Complete Cordless Obturation System

Buy a Dia-Duo, Get 1 box of Duo-Gun tips, 2 bottles of GP Obturators, 1 Duo-Pen Heating Condenser FREE!*



Duo-Pen Cordless Warm Vertical Compaction Device

- Quick heating tip reaches 230°C within 1 sec
- 5 color-coded tip sizes:
- #.04-40, .04-45, .04-50, .06-50, .08-55
- Tightly compacts and seals canals incl. lateral canals
- 3 temp. settings: 170°C, 200°C, and 230°C
- Simple operation with only two buttons
- One full charge: operation time of 500 x 5 sec.
- Replaceable Li-ion battery

*End user specials. Dealers to ship out free goods. Mix & match okay. Free goods will be the lowest priced item purchased.