



SPECIALS

APRIL | MAY | JUNE 2019

COMPOSITES | GLASS IONOMERS
TOOTH WHITENING | EQUIPMENT
CEMENTS | ADHESIVES
ETCHANTS | SEALANTS
AMALGAMS | ACCESSORIES

WHITENING

POLA

ADVANCED TOOTH
WHITENING SYSTEM

POLA LUMINATE

BRUSH ON & GO!

- No trays required - quick and easy professional brush on tooth whitener
- Apply twice a day for just 30 minutes
- 6% hydrogen peroxide, low viscosity fast drying gel
- Over 60 applications in a 5.5mL tube



POLA OFFICE+

THE WORLD'S FASTEST BLEACH

- Can be used with or without a light
- Fast and simple to use: 3 x 8 minute applications

POLA® DAY

HYDROGEN PEROXIDE BASED ADVANCE
TRAY TOOTH WHITENING SYSTEM

- From 30 minutes once a day
- Available in 3%, 7.5% and 9.5% hydrogen peroxide
- High water content
- Neutral pH and desensitizing additives

POLA® NIGHT

CARBAMIDE PEROXIDE BASED ADVANCED
TRAY TOOTH WHITENING SYSTEM

- From 45 minutes once a day
- Available in 10%, 16% and 22% carbamide peroxide
- High water content
- Neutral pH and desensitizing additives

POLA® DAY CP

SUPER FAST TRAY TOOTH WHITENING SYSTEM

- From just 15 minutes per day tray wear time
- 35% carbamide peroxide
- Contains fluoride
- Neutral pH
- Desensitizing additives

ALL POLA
WHITENING

3+1* BUY 3,
GET 1 FREE

COMPOSITES

3+1* AURA & WAVE
BUY 3, GET 1 FREE

AURA EASY

STOCK LESS, DELIVER MORE

- 4 shades cover 8 VITA classical A1-D4® shades
- Increased wear resistance and strength
- Smooth high gloss finish every time
- Non-sticky consistency makes it simple to sculpt
- Low volumetric shrinkage to offset marginal stress



AURA BULK FILL

THE SIMPLICITY OF BULK FILLING

- One bulk fill placement up to 5mm depth of cure
- 20 sec. LED light cure
- Natural-looking finish with just one universal shade
- High depth of cure without the need for a capping layer
- Easy placement in a matter of seconds



The Dental Advisor, Vol. 33,
No. 07, September 2016



WAVE, WAVE MV AND WAVE HV

FLUORIDE RELEASING, FLOWABLE
RADIOPAQUE COMPOSITES

- High compressive strength
- Nano technology
- Three viscosities
- Natural fluorescence
- No after cure shade change



ACCESSORIES

POINTS

DISPOSABLE BRUSH APPLICATORS

- Solvent resistant, non-absorbent, non-linting fibers
- Boxes of 400 in medium, fine and superfine



3+2* BUY 3,
GET 2 FREE

EQUIPMENT

RADII-CAL

LIGHTWEIGHT, RADICALLY POWERFUL LED CURING LIGHT

- Pulse technology = High intensity 1200mW/cm² plus low heat emission
- Stays charged for 1200 ten second cures (equivalent to 3h 20 min continuous cure)



2+1*
BUY 2, GET 1 FREE
RADII CAL AND
ACCESSORIES

RADII PLUS

NOT JUST AN LED CURING LIGHT

- 5 year warranty[^]
- Pulse technology = high intensity 1500mW/cm² plus low heat emission
- Stays charged for 1200 ten second cures (equivalent to 3h 20 min continuous cure)



2+1* BUY 2, GET 1 FREE
RADII PLUS AND
ACCESSORIES

[^]5 year warranty on all parts except battery and full arch bleaching LED attachment which is 2 years.

COMPOSITES | GLASS IONOMERS | TOOTH WHITENING | EQUIPMENT | CEMENTS
ADHESIVES | ETCHANTS | SEALANTS | AMALGAMS | ACCESSORIES

Redeem your free goods online using this link:
<http://www.sdi.com.au/freegoods-us>

*Please allow 2-4 weeks for delivery. Please note offers cannot be combined. Offers cannot be used in conjunction with contract pricing. Valid on invoices dated April 1st - June 30th, 2019. All claims for free goods must be made within 60 days from the end of promotional date. Free products must be of equal or lesser value. *Pota Office+ The world's fastest bleach. Based on the total treatment time compared to that of all other competitor's printed literature. Offers are valid in the US only.



YOUR
SMILE.
OUR
VISION.

SDI (North America) Inc.
1279 Hamilton Parkway
Itasca, IL 60143
Phone 1 800 228 5166
Fax: 630-361-9222
usa.canada@sdi.com.au
www.sdi.com.au